


Vascular Intervention // **Coronary**

NHancer Pro X

Microcatheter



NHancer Pro X

Microcatheter

Vascular
Intervention
Coronary



Technical Data	A-NX3	A-NX6	R-NX3	R-NX6
Guide wire size	0.014"	0.014"	0.014"	0.014"
Usable length	135 cm	135 cm	155 cm	155 cm
Hydrophilic coating length	25 cm	25 cm	60 cm	60 cm
Proximal shaft outer diameter	3.0F	2.6F	3.0F	2.6F
Distal shaft outer diameter	2.0F	2.3F	2.0F	2.3F
Proximal shaft inner diameter	0.52 mm	0.52 mm	0.52 mm	0.52 mm
Distal shaft inner diameter	0.45 mm	0.48 mm	0.45 mm	0.48 mm

Ordering Information	Configuration	Ordering number	IMDS article number	Usable length
	A-NX3	451672	NX31413525	135 cm
	A-NX6	451673	NX61413525	135 cm
	R-NX3	451674	NX31415560	155 cm
	R-NX6	451675	NX61415560	155 cm

NHancer Pro X is a trademark of IMDS.
Distributed by BIOTRONIK in selected countries.

NHANCER PRO_x

Manufacturer:
IMDS Operations B.V.
Ceintuurbaan Noord 150
9301 NZ Roden,
The Netherlands
Tel +31 (0)50 8200230
info@imds.nl

Distributor:
BIOTRONIK AG
Ackerstrasse 6
8180 Bülach, Switzerland
Tel +41 (0) 44 8645111
Fax +41 (0) 44 8645005
info.vi@biotronik.com
www.biotronik.com



© 2023 BIOTRONIK AG – All rights reserved.
Specifications are subject to modification, revision and improvement.

